


HighRanger Plus 100bp DNA Ladder (100bp - 10,000bp)

Catalogue Number: 12000

100 gel lanes





NORGEN
BIOTEK  CORPORATION

HighRanger Plus 100bp DNA Ladder (100bp - 10,000bp)

Catalogue Number: 12000

(100 gel lanes)

Description:

High quality ready-to-use double-stranded DNA ladder derived by recombinant technology for precise on-gel size and mass determination.

Features:

- Precise band sizes and masses for densitometric quantification.
- Nineteen discrete bands, ranging from 100bp to 10,000bp.
- Higher intensity reference bands at 500bp, 1,000bp and 5,000bp.

Contents:

1mL of premixed DNA ladder (0.5µg/10µL) in loading buffer (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, and 0.17% SDS).

Instructions:

Mix thoroughly. For best results, load 10µL of DNA ladder per well. For precise mass determination with a densitometer, stain gel after electrophoresis using 0.5µg/mL ethidium bromide for 30-40 minutes. The accompanying diagram shows the size and mass for each band based on 10µL ladder per well.

Storage:

Stable at room temperature. For longer term storage, -20°C is recommended.

This ladder was standardized using 10µL of DNA per lane on a 0.8 cm thick, 13 x 15 cm, 1.0% agarose gel run in TAE buffer.

Tech support: techsupport@norgenbiotek.com

Orders: orders@norgenbiotek.com

Toll Free: 1-866-667-4362

Fax: 905-227-1061

www.norgenbiotek.com