



**Urine Preservative Single Dose**  
Product # 18126

**Product Insert**

**Norgen's Urine Preservation contained in single dose ampules provides a rapid and simple method for the preservation and shipping of urine samples.**

**Intended Use**

Norgen's Urine Preservative is an aqueous preservative that allows for the rapid preservation of nucleic acids and proteins in fresh urine samples. The Urine Preservative eliminates the need to immediately process or freeze samples and allows the samples to be shipped to centralized testing facilities at ambient temperatures. Simply collect a fresh urine sample into the container of your choice, cut the end from the Urine Preservative Single Dose ampule, add the preservative to the urine sample and mix. The components of the Urine Preservative allow samples to be stored for over 2 years at room temperature with no detected degradation of urine DNA, RNA or proteins.

**Nucleic Acid Isolation From Preservative**

Urine samples collected and preserved in Norgen's Urine Preservative are compatible with most DNA and RNA isolation methods. Samples stored in the tubes have been used successfully with all of Norgen Biotek's urine-based DNA, RNA and protein isolation kits and reagents.

**Kit Components**

Component	Product # 18126
Urine Preservative Single Dose	50
ID label	50
Product Insert	1

**Advantages**

- No need to immediately process samples
- Nucleic acid and protein preservation at room temperature for over 2 years
- Ship urine samples at room temperature safely
- Compatible with most DNA/RNA/Protein isolation methods

**Specifications**

Kit Specifications	
Minimum Urine Input	5 mL
Maximum Urine Input	50 mL
Stability of DNA, RNA, Proteins at Room Temperature	Over 2 years

### **Storage Conditions and Product Stability**

All kit components can be stored for 2 years after the date of shipment without showing any reduction in performance.

### **Quality Control**

In accordance with Norgen's Quality Management System, Norgen's Urine Preservation undergoes quality assurance testing to verify that performance does not vary from lot to lot.

## **Procedure**

### **Notes Prior to Use**

- Urine of all human and animal subjects is considered potentially infectious. All necessary precautions recommended by the appropriate authorities in the country of use should be taken when working with urine.
  - Dispose of unused reagents and human specimens according to local, provincial or federal regulations.
- 
1. Collect 5 to 50 mL of midstream urine in a sterile container of your choice. First flow urine is not recommended for collection.
  2. Carefully cut the tip off the Urine Preservative Single Dose ampule and dispense the contents into the container with the urine sample. Tighten container lid.
  3. Invert gently several times to mix preservative with urine sample.
  4. Fill in any required information on the provided ID Labels and affix to the container with the preserved urine sample.
  5. Store or ship the urine at room temperature.
  6. RNA, DNA and proteins can be purified from the preserved urine sample using a number of different methods, including Norgen's urine-based sample preparation kits.

### **Technical Support**

Contact our Technical Support Team between the hours of 9:00 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362.

Technical support can also be obtained from our website ([www.norgenbiotek.com](http://www.norgenbiotek.com)) or through email at [support@norgenbiotek.com](mailto:support@norgenbiotek.com).

3430 Schmon Parkway, Thorold, ON Canada L2V 4Y6  
Phone: (905) 227-8848  
Fax: (905) 227-1061  
Toll Free in North America: 1-866-667-4362