

3430 Schmon Parkway Thorold, ON, Canada L2V 4Y6 Phone: 866-667-4362 • (905) 227-8848 Fax: (905) 227-1061

Email: techsupport@norgenbiotek.com

Total RNA Purification Kit (Magnetic Bead System) Product Insert Product # 75400

Norgen's Total RNA Purification Kit (Magnetic Bead System) offers a fast and reliable method for isolating and purifying total RNA from a wide range of sample types, including cultured animal cells, tissue samples, blood, plasma, serum, bacteria, yeast, fungi, plants, and viruses. The kit efficiently purifies RNA of all sizes, from large mRNA and ribosomal RNA to microRNA (miRNA) and small interfering RNA (siRNA). The purification process selectively isolates RNA from other cellular components, such as proteins, without requiring phenol or chloroform. The resulting RNA is of exceptional integrity and is suitable for various downstream applications, including real-time PCR, reverse transcription PCR, next-generation sequencing (NGS), Northern blotting, RNase protection assays, primer extension, and expression array analyses.

Norgen's Purification Technology

Purification relies on magnetic beads that selectively bind RNA under optimized conditions. The process begins with lysing cells or tissues using Buffer RL. Isopropanol is added to the lysate, enabling RNA to bind to the magnetic beads based on ionic concentrations. The binding chemistry ensures only RNA attaches to the beads, while contaminating proteins are removed in the supernatant. The bound RNA is then washed with Wash Solution A to eliminate residual impurities, followed by elution of the purified RNA with Elution Solution A. The resulting RNA is highly pure and suitable for various downstream applications.

Kit Components

Component	Product #75400 (50 samples)
Buffer RL	70 mL
Magnetic Beads C	2.4 mL
Wash Solution A	38 mL
Elution Solution A	6 mL
Elution Tubes (1.7 mL)	50
Product Insert	1

Advantages

- Fast, reproducible and easy processing of samples.
- Isolates total RNA including microRNAs.
- Isolates RNA from a variety of microbial samples (Bacteria, Viruses, Yeast).
- Isolates RNA from cells and a variety of tissues (Lung, Heart, Spleen).
- Tested on a variety of bodily fluids (Blood, Plasma, Saliva, Swabs).
- Magnetic Bead protocol can easily be adapted on a variety of automation systems (IsoPure[™], KingFisher, Hamilton Star/Vantage, and Tecan)
- Purified RNA is suitable for a variety of downstream applications, including qPCR or NGS sequencing.

Specifications

Kit Specifications			
Number of Preps	50		
Size of Purified RNA	All sizes, including Small RNA (<200 nt)		
Maximum Amount of Starting Material:			
Animal Cells	3 x 10 ⁶ cells		
Animal Tissues	10 mg (for most tissues*)		
Blood	100 μL		
Plasma/Serum	200 μL		
Bacteria	1 x 10 ⁹ cells		
Yeast	1 x 10 ⁸ cells		
Fungi	50 mg		
Plant Tissues	50 mg or ≤5x10 ⁶ plant cells		
Virus	100 μL viral suspension		
Average Yield *	4 – 20 μg		
Average Purity (OD260/280)	1.8 – 2.0		
Time to Complete 50 Purifications	15 minutes (hands on time)		
(automated)	15 minutes (hands-on time)		
Time to Complete 50 Purifications	40 minutes (hands-on time)		
(manual)	40 minutes (nanus-on time)		

^{*} Average RNA yield will vary depending on the sample type

Storage Conditions and Product Stability

All solutions should be kept tightly sealed and stored at room temperature. These kits are stable for 2 years after the date of shipment.

Precautions and Disclaimers

This kit is designed for research purposes only. It is not intended for human or diagnostic use. Ensure that a suitable lab coat, disposable gloves and protective goggles are worn when working with chemicals. For more information, please consult the appropriate Safety Data Sheets (SDS). These are available as convenient PDF files online at www.norgenbiotek.com.

The **Buffer RL** contains guanidinium salts and should be handled with care. Guanidinium salts form highly reactive compounds when combined with bleach, thus care must be taken to properly dispose of any of these solutions.

Blood of all human and animal subjects is considered potentially infectious. All necessary precautions recommended by the appropriate authorities in the country of use should be taken when working with whole blood.

Customer-Supplied Reagents and Equipment

For All Protocols

- Magnetic Bead Separation Stand (For manual Isolation)
- IsoPureTM 96 or IsoPureTM Mini Purification System (For automated Isolation)
- Micropipettors

- Microcentrifuge tubes
- β-mercaptoethanol (optional)
- Isopropanol
- Temperature adjustable (50°C) incubator
- Nuclease-free water
- DNase I (optional)
- Phosphate Buffered Saline (PBS) (for Cultured Animal Cells protocol)
- Cell Disruption Tools such as mortar and pestle, rotor-stator homogenizer (for Animal Tissue protocol)
- Syringe with a 22G needle (for Animal Tissue protocol)
- Tube Revolver Rotator

For Animal Cell Protocol

PBS (RNase-free)

For Animal Tissue Protocol

- Liquid nitrogen
- Mortar and pestle

For Nasal or Throat Swabs

• Sterile, single-use cotton swabs

For Bacterial Protocol

- Lysozyme-containing TE Buffer:
- For Gram-negative bacteria, 1 mg/mL lysozyme in TE Buffer
- For Gram-positive bacteria, 3 mg/mL lysozyme in TE Buffer

For Yeast Protocol

- Resuspension Buffer with Lyticase:
- 50 mM Tris pH 7.5
- 10 mM EDTA
- 1 M Sorbital
- 1 unit/L Lyticase

For Fungi Protocol

- Liquid nitrogen
- Mortar and pestle

For Plant Protocol


- Liquid nitrogen
- Mortar and pestle

For Plasma/Serum Protocol

MS2 RNA (0.8 μg/μl). (Roche, Cat. No. 10165948001)

Flow Chart

Procedure for using Norgen's Total RNA Purification Kit (Magnetic Bead System) –IsoPure™ 96




Run Instrument




Flow Chart

Procedure for using Norgen's Total RNA Purification Kit (Magnetic Bead System) −IsoPure™ Mini



Flow Chart

Procedure for Norgen's Total RNA Purification Kit (Magnetic Bead System) - Manual Method



Procedures

Working with RNA

RNases are very stable and robust enzymes that degrade RNA. Autoclaving solutions and glassware are not always sufficient to actively remove these enzymes. The first step when preparing to work with RNA is to create an RNase-free environment. The following precautions are recommended as your best defense against these enzymes.

- The RNA area should be located away from microbiological workstations.
- Clean, disposable gloves should be worn at all times when handling reagents, samples, pipettes, disposable tubes, etc. It is recommended that gloves are changed frequently to avoid contamination.
- There should be designated solutions, tips, tubes, lab coats, pipettes, etc. for RNA only.
- All RNA solutions should be prepared using at least 0.05% DEPC-treated autoclaved water or molecular biology grade nuclease-free water.
- Clean all surfaces with commercially available RNase decontamination solutions.
- When working with purified RNA samples, ensure that they remain on ice during downstream applications.

- The procedure provides information and steps to follow for extraction of Total RNA using an Automated as well as a Manual method.
- In each method, the steps for preparing the lysate are different depending on the starting material (Section 1). However, the subsequent steps from all lysate types for automation as well as manual extraction are mentioned in Section 2.
- Please ensure that the correct procedure for preparing the lysate from your starting material is followed.
- Ensure that all solutions are at room temperature prior to use, and that no precipitates have formed. If necessary, warm the solutions and mix well until the solutions become clear again.
- Prepare a working concentration of the Wash Solution A by adding:
 - a) 90 mL of 96 -100% ethanol (provided by the user) to each of the bottles containing 38 mL of concentrated **Wash Solution A**. This will give a final volume of 128 mL.
- The labels on the bottles have a box that may be checked to indicate that the ethanol has been added.
- Always vortex the RNA Magnetic Bead C before use.
- Optional: The use of β-mercaptoethanol in lysis is highly recommended for most animal tissues, particularly those known to have a high RNAse content (ex: pancreas), as well as for most plant tissues and nasal and throat swabs. It is also recommended for users who wish to isolate RNA for sensitive downstream applications. Add 10 μL of β-mercaptoethanol (provided by the user) to each 1 mL of Buffer RL required. β-mercaptoethanol is toxic and should be dispensed in a fume hood. Alternatively, the Buffer RL can be used as provided.
- It is important to work quickly during this procedure.
- Preheat the incubator(s) according to the temperatures required (37°C or 55°C or 65°C). For automation procedures, the automation platform that is being used may have a heat block to provide appropriate temperatures for incubation and thus an incubator might not be required.

Total RNA Isolation using Automation and Manual Method

This section outlines step-by-step instructions for extracting Total RNA from various sample types using both automated and manual methods. The automation method provides guidance for using IsoPure™ instruments, but the protocols can also be easily adapted for other magnetic bead-based automation platforms, such as KingFisher, Hamilton, or Tecan.

Section 1 details the preparation of lysates from various sample types, applicable to both automated and manual RNA extraction methods.

Section 2 is divided into multiple subsections, offering detailed instructions for extracting RNA from the lysates prepared in **Section 1**, using either automated systems or manual procedures.

To prepare for Total RNA extraction using automated systems, please follow the steps outlined below.

Using the IsoPure™ Mini or IsoPure™ 96 Automation System

The procedure outlined below is for extraction of Total RNA using Automation platforms IsoPure ™ 96 (Cat. AP1096) and IsoPure ™ Mini (Cat. AP1016).

Set-up the IsoPure[™] Plates

Notes Prior to Use

- The numbering of plates is based on where it should be placed on the carousel in the instrument.
- IsoPure[™] 96 uses multiple plates for RNA extraction while IsoPure[™] Mini uses a single plate for RNA extraction.
- For IsoPure[™] 96 most of the lysate preparation is performed in plate 2 while other plates are used for performing washing and elution steps.
- For IsoPure™ Mini the lysate preparation is performed in column 1 and/or 7 while washing. Elutions steps could be performed in other columns of the same plate.
- For IsoPure[™] 96, RNA elution step is always performed in plate 8.
- Prepare the plates and label them based on the sample type as mentioned in the tables below.

Table 1. Set-up for Total RNA extraction from All Sample Types

Prepare plates as per the volumes mentioned in the table below. Detailed description of step-wise process can be found in **Section 1 and 2**.

Component	IsoPure™ Mini	IsoPure™ 96	Location on Carousel	Amount
Tip Comb	In the instrument	Plate 1	1	Empty Plate
Lysate Column 1 / Column 7		Plate 2	2	Variable
Wash Solution A Column 3 / Column 9,		Plate 4, Plate 5, &	4, 5 & 6	600 µL
	Column 4 / Column 10 and	Plate 6		
	Column 5 / Column 11			
Elution Solution A Column 6 / Column 12		Elution Plate	8	50 – 100 μL

Norgen's Total RNA Purification Kit (Magnetic Bead System) isolates total RNA with minimal amounts of genomic DNA contamination. However, an optional **DNA Removal Protocol** is provided for maximum removal of residual DNA that may affect sensitive downstream applications. It is recommended that

Norgen's RNase-Free DNase I Kit (Product # 25710) be used for this step. Prepare a DNase working solution as mentioned below when using DNase treatment with the extraction.

i. For every DNase reaction to be performed, prepare a DNase I Mix of 15 μL of DNase I and 100 μL of Enzyme Incubation Buffer A using Norgen's RNase-Free DNase I Kit (Product# 25710). Mix gently by inverting the tube a few times. DO NOT VORTEX.

Note: If using an alternative DNase I, prepare a working stock of 0.25 Kunitz/ μ L RNase-free DNase I solution according to the manufacturer's instructions. A 100 μ L aliquot is required for each sample to be treated. Set-up the plates as per **Table 2** when performing DNase treatment of the samples.

Table 2. Set-up for Total RNA extraction from All Sample Types – with DNase Treatment

Prepare plates as per the volumes mentioned in the table below. Detailed description of step-wise process can be found in **Section 1 and 2.**

Component	IsoPure™ Mini	IsoPure™ 96	Location on Carousel	Amount
Tip Comb	In the instrument	Plate 1	1	Empty Plate
Lysate	Lysate Column 1 / Column 7		2	Variable
DNase I Mix Column 2 / Column 8		Plate 3	3	100 μL
Wash Solution A Column 3 / Column 9,		Plate 4, Plate 5, &	4, 5 & 6	600 µL
	Column 4 / Column 10 and	Plate 6		
	Column 5 / Column 11			
Elution Solution A Column 6 / Column 12		Elution Plate	8	50 – 100 μL

Section 1. Lysate Preparation from Various Sample Types

Notes Prior to Use

- Fresh or frozen tissues may be used for the procedure.
- Tissues should be immediately frozen and stored at -20°C or -70°C. Tissues may be stored at -70°C for several months.
- It is recommended that no more than 20 mg of tissue be used.

1A. Lysate Preparation from Cultured Animal Cells

- Cells grown in suspension or monolayer may be used.
- The maximum recommended input of cells is 3 x 10⁶. A hemocytometer can be used in conjunction with a microscope to count the number of cells. As a general guideline, a confluent 3.5 cm plate of HeLa cells will contain 10⁶ cells.
- Cell pellets can be stored at -70°C for later use or used directly in the procedure. Determine the number of cells present before freezing.
- Frozen cell pellets should not be subjected to repeated cycles of freeze and thaw prior to beginning the protocol.
- Frozen cell pellets should not be thawed prior to beginning the protocol. Add the Buffer RL directly to the frozen cell pellet (Step 1A (ii) c).

1A (i). Cell Lysate Preparation from Cells Growing in a Monolayer

- a) Aspirate media and wash cell monolayer with an appropriate amount of PBS. Aspirate PBS.
- b) Add 350 µL of **Buffer RL** directly to culture plate.
- c) Lyse cells by gently tapping culture dish and swirling buffer around plate surface for five minutes.
- d) Transfer lysate to **Plate 2** (IsoPure[™] 96), **Column 1 and/or Column 7** (IsoPure[™] Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction.
- e) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

Note: For input amounts greater than 10⁶ cells, it is recommended that the lysate is passed through a 25-gauge needle attached to a syringe 5-10 times at this point, in order to shear the genomic DNA prior to loading into the well or column.

1A (ii). Cell Lysate Preparation from Cells Growing in Suspension and Lifted Cells

- a) Transfer cell suspension to an RNase-free tube (not provided) and centrifuge at no more than 200 x g (~2,000 RPM) for 10 minutes to pellet cells.
- b) Carefully decant the supernatant. A few μ L of media may be left behind with the pellet in order to ensure that the pellet is not dislodged.
- c) Add 350 µL of **Buffer RL** to the pellet and dissolve the pellet by mixing with a pipette. Ensure that the entire pellet is completely dissolved before proceeding to the next step.
- d) Transfer lysate to **Plate 2** (IsoPure[™] 96), **Column 1 and/or Column 7** (IsoPure[™] Mini) or RNase-free **microcentrifuge tube** (not provided) for manual extraction.
- e) Proceed to **Section 2:** RNA Isolation from all Types of Lysates.

Note: For input amounts greater than 10⁶ cells, it is recommended that the lysate is passed through a 25-gauge needle attached to a syringe 5-10 times at this point, in order to shear the sample contents.

If the volume of cell suspension is 1 mL or less, **steps a − c** can be performed in **Plate 2** (IsoPureTM 96), **Column 1 and/or Column 7** (IsoPureTM Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction. Make sure to seal the plate with an Adhesive Film (not provided). In this case, proceed to **Section 2** directly after performing **step c**.

1B. Lysate Preparation from Animal Tissues

- Norgen's Total RNA Purification Kit (Magnetic Bead System) is designed for isolating RNA from small amounts of tissue sample (up to 10 mg in most cases). If a larger amount of starting material is desired, Norgen's Animal Tissue RNA Purification Kit (Cat. #25700) should be used.
- RNA in animal tissues is not protected after harvesting until it is disrupted and homogenized.
 Thus, it is important that the procedure is carried out as quickly as possible, particularly the Cell Lysate Preparation step.
- Fresh or frozen tissues may be used for the procedure. Tissues should be flash-frozen in liquid nitrogen and transferred immediately to a -70°C freezer for long-term storage. Tissues may be stored at -70°C for several months. When isolating total RNA from frozen tissues ensure that the tissue does not thaw during weighing or prior to grinding with the mortar and pestle.

- Tissues stored in RNA stabilization reagents such as RNAlater® are compatible with this isolation procedure. Prior to isolation, carefully remove the tissue from the storage reagent using forceps, and dry excessive liquid.
- The maximum recommended input of tissue varies depending on the type of tissue being used.
 Please refer to **Table 3** below as a guideline for maximum tissue input amounts. If your tissue of interest is not included in the table below, we recommend starting with an input of no more than 10 mg.

Table 3. Recommended Maximum Input Amounts of Different Tissues

Tissue	Maximum Input Amount
Brain	25 mg
Heart	5 mg
Kidney	10 mg
Liver	10 mg
Lung	10 mg
Spleen	10 mg

1B. Cell Lysate Preparation from Animal Tissues

- a) Excise the tissue sample from the animal.
- b) Determine the amount of tissue by weighing. Please refer to **Table 3** for the recommended maximum input amounts of different tissues. For tissues not included in the table, we recommend starting with an input of no more than 10 mg.
- c) Transfer the tissue into a mortar that contains an appropriate amount of liquid nitrogen to cover the sample. Grind the tissue thoroughly using a pestle.
- d) Allow the liquid nitrogen to evaporate, without allowing the tissue to thaw.
- e) Add 600 μL of **Buffer RL** to the tissue sample and continue to grind until the sample has been homogenized. Homogenize by passing the lysate 5 10 times through a 25-gauge needle attached to a syringe.
- f) Using a pipette, transfer the lysate into an RNase-free microcentrifuge tube (not provided).
- g) Spin the lysate for 2 minutes to pellet any cell debris and transfer 500 μL lysate to **Plate 2** (IsoPureTM 96), **Column 1 and/or Column 7** (IsoPureTM Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction.
- h) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

1C. Lysate Preparation from Blood

- This procedure is for the isolation of RNA from whole blood. For the isolation of RNA from plasma or serum samples, please see Section 1J.
- Blood of all human and animal subjects is considered potentially infectious. All necessary
 precautions recommended by the appropriate authorities in the country of use should be taken
 when working with whole blood.
- It is recommended that no more than 100 μL of blood be used in order not to saturate the magnetic beads.
- We recommend the use of this kit to isolate RNA from non-coagulating fresh blood using EDTA as the anti-coagulant.
- It is important to work quickly during this procedure.

1C. Cell Lysate Preparation from Blood

- a) Transfer up to 100 µL of non-coagulated blood to **Plate 2** (IsoPure[™] 96), **Column 1 and/or Column 7** (IsoPure[™] Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction.
- b) Add 350 μL of Buffer RL to the blood.
- c) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

1D. Lysate Preparation from Nasal or Throat Swabs

Notes Prior to Use

- Body fluid of all human and animal subjects is considered potentially infectious. All necessary
 precautions recommended by the appropriate authorities in the country of use should be taken
 when working with these samples.
- It is important to work quickly during this procedure.

1D. Lysate Preparation from Nasal or Throat Swabs

- a) Add 600 µL of Buffer RL to an RNase-free microcentrifuge tube (not provided).
- b) Gently brush a sterile, single-use cotton swab inside the nose or mouth of the subject.
- c) Using sterile techniques, cut the cotton tip where the nasal or throat cells were collected and place into the microcentrifuge tube containing the **Buffer RL**. Close the tube. Vortex gently and incubate for 5 minutes at room temperature.
- i) Transfer 500 μL lysate to **Plate 2** (IsoPureTM 96), **Column 1 and/or Column 7** (IsoPureTM Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction.

Note: In case of swab samples collected in Norgen's Total Nucleic Acid Preservation Tubes (Cat# 69200) transfer 500 μL preserved sample directly to Plate 2 (IsoPureTM 96), Column 1 and/or Column 7 (IsoPureTM Mini) or the 1.5 mL microcentrifuge tube (Manual extraction) and proceed to Section 2.

j) Proceed to **Section 2**: RNA Isolation from All Types of Lysates.

1E. Lysate Preparation from Bacteria

- Prepare the appropriate lysozyme-containing TE Buffer as indicated in **Table 4.** This solution should be prepared with sterile, RNase-free TE Buffer, and kept on ice until needed. These reagents are to be provided by the user.
- It is recommended that no more than 10⁹ bacterial cells be used in this procedure. Bacterial growth can be measured using a spectrophotometer. As a general rule, an *E. coli* culture containing 1 x 10⁹ cells/mL has an OD₆₀₀ of 1.0.
- For RNA isolation, bacteria should be harvested in log-phase growth.
- Bacterial pellets can be stored at -70°C for later use, or used directly in this procedure.
- Frozen bacterial pellets should not be thawed prior to beginning the protocol. Add the Lysozymecontaining TE Buffer directly to the frozen bacterial pellet (Step 1E. c).

1E. Cell Lysate Preparation from Bacteria

- a) Pellet bacteria by centrifuging at 14,000 x g (~14,000 RPM) for 1 minute. (*Note:* Centrifuge at 4000 RPM for 3 minutes while centrifuging bacterial suspension in a 96-well plate).
- b) Decant supernatant, and carefully remove any remaining media by aspiration.
- c) Resuspend the bacteria thoroughly in 100 µL of the appropriate lysozyme-containing TE buffer by vortexing. Incubate at room temperature for the time indicated in **Table 4**.
- d) Add 350 μ L of **Buffer RL** to the pellet and dissolve the pellet by mixing with a pipette. Ensure that the entire pellet is completely dissolved before proceeding to the next step.
- e) Transfer lysate to **Plate 2** (IsoPure[™] 96), **Column 1 and/or Column 7** (IsoPure[™] Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction.

Note: If the volume of bacterial suspension is 1 mL or less, steps a − d can be performed in Plate 2 (IsoPureTM 96), Column 1 and/or Column 7 (IsoPureTM Mini) or the RNase-free microcentrifuge tube (not provided) for manual extraction. Make sure to seal the plate with an Adhesive Film (not provided). In this case proceed to Section 2 directly after performing step d.

f) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

Table 4: Incubation Time for Different Bacterial Strains

Bacteria Type	Lysozyme Concentration in TE Buffer	Incubation Time
Gram Negative	1 mg/mL	5 min
Gram Positive	3 mg/mL	10 min

1F. Lysate Preparation from Yeast

Notes Prior to Use

- Prepare the appropriate amount of Lyticase-containing Resuspension Buffer, considering that 100 μ L of buffer is required for each preparation. The Resuspension Buffer should have the following composition: 50 mM Tris, pH 7.5, 10 mM EDTA, 1M Sorbital, 0.1% β -mercaptoethanol and 1 unit/ μ L Lyticase. This solution should be prepared with sterile, RNAse-free reagents, and kept on ice until needed. These reagents are to be provided by the user.
- It is recommended that no more than 10⁷ yeast cells or 1 mL of culture be used for this procedure.
- For RNA isolation, yeast should be harvested in log-phase growth.
- Yeast can be stored at -70°C for later use, or used directly in this procedure.
- Frozen yeast pellets should not be thawed prior to beginning the protocol. Add the Lyticase-containing Resuspension Buffer directly to the frozen yeast pellet (Step 1F. c).

1F. Cell Lysate Preparation from Yeast

- a) Pellet yeast by centrifuging at 14,000 x g (~14,000 RPM) for 1 minute. (*Note:* Centrifuge at 4000 RPM for 3 minutes while centrifuging bacterial suspension in a 96-well plate).
- b) Decant supernatant, and carefully remove any remaining media by aspiration.
- c) Resuspend the yeast thoroughly in 100 μL of Lyticase-containing Resuspension Buffer by vortexing.
- d) Incubate at 37°C for 10 minutes.

- e) Add 350 µL of **Buffer RL** to the pellet and dissolve the pellet by mixing with a pipette. Ensure that the entire pellet is completely dissolved before proceeding to the next step.
- f) Transfer lysate to **Plate 2** (IsoPure[™] 96), **Column 1 and/or Column 7** (IsoPure[™] Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction.

Note: If the volume of yeast suspension is 1 mL or less, step a − e can be performed in Plate 2 (IsoPureTM 96), Column 1 and/or Column 7 (IsoPureTM Mini) or the RNase-free microcentrifuge tube (not provided) for manual extraction. Make sure to seal the plate with an Adhesive Film (not provided). In this case proceed to Section 2 directly after performing step e.

g) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

1G. Lysate Preparation from Fungi

Notes Prior to Use

- Fresh or frozen fungi may be used for this procedure. Fungal tissues should be flash frozen in liquid nitrogen and transferred immediately to a -70°C freezer for long-term storage. Fungi may be stored at -70°C for several months. Do not allow frozen tissues to thaw prior to grinding with the mortar and pestle in order to ensure that the integrity of the RNA is not compromised.
- It is recommended that no more than 50 mg of fungi be used for this procedure in order not to saturate the magnetic beads.

1G. Cell Lysate Preparation from Fungi

- a) Determine the amount of fungi by weighing. It is recommended that no more than 50 mg of fungi be used for the protocol.
- b) Transfer the fungus into a mortar that contains an appropriate amount of liquid nitrogen to cover the sample. Grind the fungus thoroughly using a pestle.

Note: At this stage the ground fungus may be stored at -70°C, such that the RNA purification can be performed at a later time.

- c) Allow the liquid nitrogen to evaporate, without allowing the tissue to thaw.
- d) Add 600 µL of **Buffer RL** to the tissue sample and continue to grind until the sample has been homogenized.
- e) Using a pipette, transfer the lysate into an RNase-free microcentrifuge tube (not provided).
- f) Spin the lysate for 2 minutes to pellet any cell debris. Transfer 500 μL lysate to **Plate 2** (IsoPureTM 96), **Column 1 and/or Column 7** (IsoPureTM Mini) or the RNase-free **microcentrifuge tube** (not provided) for manual extraction.
- g) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

1H. Lysate Preparation from Plant

Notes Prior to Use

• The maximum recommended input of plant tissue is 50 mg or 5 x 10⁶ plant cells.

- Both fresh and frozen plant samples can be used for this protocol. Samples should be flash-frozen in liquid nitrogen and transferred immediately to a -70°C freezer for long-term storage. Do not allow frozen tissues to thaw prior to grinding with the mortar and pestle in order to ensure that the integrity of the RNA is not compromised.
- It is important to work quickly during this procedure.

1H. Cell Lysate Preparation from Plant

a) Determine the amount of plant by weighing and transfer ≤50 mg plant tissue or 5 x 10⁶ plant cells into a mortar that contains an appropriate amount of liquid nitrogen to cover the sample. Grind the sample into a fine powder using a pestle in liquid nitrogen.

Note: If stored frozen samples are used, do not allow the samples to thaw before transferring to the liquid nitrogen.

- b) Allow the liquid nitrogen to evaporate, without allowing the tissue to thaw.
- c) Add 600 µL of **Buffer RL** to the tissue sample and continue to grind until the sample has been homogenized.
- d) Using a pipette, transfer the lysate into an RNase-free microcentrifuge tube (not provided). Spin the lysate for 2 minutes to pellet any cell debris. Transfer up to 500 μL supernatant to Plate 2 (IsoPureTM 96), Column 1 and/or Column 7 (IsoPureTM Mini) or the RNase-free microcentrifuge tube (not provided) for manual extraction.
- e) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

11. Lysate Preparation from Viral Suspension

Notes Prior to Use

- For the isolation of integrated viral RNA, follow Section 1A if the starting material is tissue, follow Section 1B if the starting material is cell culture, follow Section 1C if the starting material is blood, or follow Section 1D if the starting material is a nasal or throat swab.
- For the isolation of RNA from free viral particles, follow the procedure below.
- It is recommended that no more than 100 μL of viral suspension be used in order to perform effective isolation.
- It is important to work quickly during this procedure.

11. Cell Lysate Preparation from Viral Suspension

- a) Transfer 100 μL of viral suspension to **Plate 2** (IsoPureTM 96), **Column 1 and/or Column 7** (IsoPureTM Mini) or the **1.5 mL microcentrifuge tube** (Manual extraction).
- b) Add 350 µL of **Buffer RL**.
- c) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

1J. Lysate Preparation from Plasma or Serum

- Plasma or Serum of all human and animal subjects is considered potentially infectious. All
 necessary precautions recommended by the appropriate authorities in the country of use should
 be taken when working with plasma or serum.
- We recommend the use of this kit to isolate RNA from plasma or serum prepared by standard protocol from non-coagulating fresh blood using EDTA or sodium citrate as the anti-coagulant.

- Plasma prepared from fresh blood using heparin as an anti-coagulant is not suitable for use with this protocol. For heparin-prepared samples follow the protocol in section 1C, Lysate Preparation from Blood.
- It is recommended that no more than 200 µL of plasma or serum be used in order to prevent saturation of magnetic bead.
- Avoid multiple freeze-thaw cycle of the plasma or serum sample. Aliquot to the appropriate volume for usage prior to freezing.
- It is important to work quickly during this procedure.

1J. Lysate Preparation from Plasma or Serum

- a) Transfer 100 µL of plasma or serum to **Plate 2** (IsoPure[™] 96), **Column 1 and/or Column 7** (IsoPure[™] Mini) or the **1.5 mL microcentrifuge tube** (Manual extraction).
- b) Add 350 µL of Buffer RL.
- c) Optional: Add 0.7 µL of 0.8 µg/µL MS2 RNA per sample.

Note: The use of MS2 RNA could increase the consistency of downstream applications such as RT-PCR. However, the use of MS2 RNA is not recommended for applications involving global gene expression analysis such as microarrays or sequencing.

d) Proceed to **Section 2:** RNA Isolation from All Types of Lysates.

Section 2: Total RNA Purification from All Types of Lysates

2.1 Total RNA Purification using IsoPure™ 96 system

- a) Setup plates according to **Table 1**. Add 600 µL of **Wash Solution A** to **Plate 4**, **5 and 6** and keep it at **position 4**, **5 and 6** respectively on the carousel.
- b) Add $50 100 \mu L$ of **Elution Solution A** to the **Elution plate** and place it at **position 8** on the carousel.
- c) Place sample plate (Plate 2) in the machine at position 2.
- d) Place a clean tip comb in an empty plate (Plate 1) and keep it at position 1.
- e) After setting up, run the program TOTRNA96.
- f) Prepare a mixture of **Isopropanol** and 40 μL **Magnetic Beads C** as per **Table 5** for the number of samples extracted. (For example, if RNA is to be extracted from 10 samples, and if the lysate volume obtained in **Section 1** is 500 μL, mix 5 mL of **Isopropanol** with 400 μL **Magnetic Beads C** and add 540 μL of the mixture to each well of **Plate 2**). A fresh bulk mixture can be made every time RNA isolation is performed. Make 10% extra volume to avoid volume loss due to pipetting.
- g) After a few minutes of incubation, the instrument will pause and allow the user to add mixture prepared in **step f**. Open the cover and remove the plate from the instrument and add appropriate volume of **Isopropanol** (not provided) and **Magnetic Beads C** mixture that was prepared in **step f**.
- h) Put the plate back into the instrument and press OK. The instrument will continue to perform the extraction.
- i) The instrument will prompt after the procedure is complete. Remove the **Elution plate** from the instrument and the RNA is now ready for further downstream processing.
- j) The purified RNA sample may be stored at -20°C for a few days. It is recommended that samples be placed at -70°C for long-term storage.

2.2 Total RNA Purification using IsoPure™ Mini system

- a) Setup plates according to **Table 1.** Add 600 µL of **Wash Solution A** to **Column 3 and/or Column 9, Column 4 and/or Column 10, Column 5 and/or Column 11.**
- b) Add $50 100 \,\mu\text{L}$ of Elution Solution A to Column 6 and/or 12.
- c) Place a clean tip comb into IsoPure™ Mini.
- d) After setting up, place the plate into the instrument and run the program TOTRNA16.
- e) Prepare a mixture of **Isopropanol** and 40 μL **Magnetic Beads C** as per **Table 5** for the number of samples extracted. (For example, if RNA is to be extracted from 10 samples, and if the lysate volume obtained in **Section 1** is 500 μL, mix 5 mL of **Isopropanol** with 400 μL **Magnetic Beads C** and add 540 μL of the mixture to each well of **Colum 1 and/or Column 7**). A fresh bulk mixture can be made every time RNA isolation is performed. Make 10% extra volume to avoid volume loss due to pipetting.
- f) After few minutes of incubation, the instrument will pause and allow the user to add mixture prepared in step e. Open the cover and remove the plate from the instrument and add appropriate volume of Isopropanol (not provided) and Magnetic Beads C mixture that was prepared in step e to 1 Column 3 and/or Column 7.
- g) Put the plate back into the instrument and press OK. The instrument will continue to perform the extraction.
- h) The instrument will prompt after the procedure is complete. Remove the plate from the instrument and the RNA is now ready for further downstream processing.
- i) Transfer the elution from **Column 6 and/or Column 12** and transfer it to a fresh plate for storage. The purified RNA sample may be stored at -20°C for a few days. It is recommended that samples be placed at -70°C for long-term storage.

Table 5. Isopropanol and Magnetic Beads C Volume for Each Sample Type

Section	Sample Type	Lysate Volume	Isopropanol	Magnetic Beads C
Section	Sample Type	(µL)	(µL)	(μL)
1A (i)	Cells (Monolayer)	350	350	40
1A (ii)	Cells (Suspension)	350	350	40
1B	Animal Tissues	500	500	40
1C	Blood	450	450	40
1D	Nasal or Throat Swabs	500	500	40
1E	Bacteria	350	350	40
1F	Yeast	350	350	40
1G	Fungi	500	500	40
1H	Plant	500	500	40
11	Viral Suspension	450	450	40
1J	Plasma or Serum	450	450	40

2.3 Total RNA Purification using IsoPure™ 96 system – with DNase Treatment

- a) Setup plates according to Table 2. Add 100 μL of DNase I Mix to Plate 3 and place it at position 3.
- b) Add 600 μ L of **Wash Solution A** to **Plate 4, 5 and 6** and place it at **position 4, 5 and 6** respectively on the carousel.
- c) Add $50 100 \mu L$ of **Elution Solution A** to the **Elution plate** and place it at **position 8** on the carousel.

- d) Place sample plate (Plate 2) in the machine at position 2.
- e) Place a clean tip comb in an empty plate (Plate 1) and keep it at position 1.
- f) After setting up, run the program TOTRNA96DNA.
- g) Prepare a mixture of **Isopropanol** (not provided) and 40 μL of **Magnetic Beads C** as per **Table 5** for the number of samples extracted. (For example, if RNA is to be extracted from 10 samples, and if the lysate volume obtained in **Section 1** is 500 μL, mix 5 mL of **Isopropanol** (not provided) with 400 μL **Magnetic Beads C** and add 540 μL of the mixture to each well of **Plate 2**). A fresh bulk mixture can be made every time RNA isolation is performed. Make 10% extra volume to avoid volume loss due to pipetting.
- h) After a few minutes of incubation, the instrument will pause and allow the user to add mixture prepared in **step g**. Open the cover and remove the plate from the instrument and add appropriate volume of **Isopropanol** (not provided) and **Magnetic Beads C** mixture that was prepared in **step g**.
- i) Put the plate back into the instrument and press OK. The instrument will continue to perform the extraction.
- j) To rebind the sample to the magnetic beads, prepare a mixture of 200 μL **Buffer RL** and 300 μL Isopropanol (not provided) for the number of samples extracted. (For example, if RNA is to be extracted from 10 samples, and mix 2 mL of **Buffer RL** with 3 mL **Isopropanol**). A fresh bulk mixture can be made every time RNA isolation is performed. Make 10% extra volume to avoid volume loss due to pipetting.
- k) After a few minutes, the instrument will pause and allow the user to add mixture prepared in **step** j. Open the cover and remove the plate from the instrument and 500 μL of **Buffer RL** and **Isopropanol** mixture that was prepared in **step** j to each well of **Plate 3**.
- I) Put the plate back into the instrument and press OK. The instrument will continue to perform the extraction.
- m) The instrument will prompt after the procedure is complete. Remove the **Elution plate** from the instrument and the RNA is now ready for further downstream processing.
- n) The purified RNA sample may be stored at -20°C for a few days. It is recommended that samples be placed at -70°C for long-term storage.

2.4 Total RNA Purification using IsoPure[™] Mini system – with DNase Treatment

- a) As per the sample being isolated, setup plates according to **Table 2.** Add 100 μ L of **DNase I Mix** to **Column 2 and/or Column 8.**
- b) Add 600 µL of Wash Solution A to Column 3 and/or Column 9, Column 4 and/or Column 10, Column 5 and/or Column 11.
- c) Add $50 100 \,\mu\text{L}$ of Elution Solution A to Column 6 and/or Column 12.
- d) Place a clean tip comb into IsoPure™ Mini.
- e) After setting up, place the plate into the instrument and run the program TOTRNA16D.
- f) Prepare a mixture of **Isopropanol** (not provided) and 40 μL **Magnetic Beads C** as per **Table 5** for the number of samples extracted. (For example, if RNA is to be extracted from 10 samples, and if the lysate volume obtained in **Section 1** is 500 μL, mix 5 mL of **Isopropanol** with 400 μL **Magnetic Beads C** and add 540 μL of the mixture to each well of **Colum 1 and/or Column 7**). A fresh bulk mixture can be made every time RNA isolation is performed. Make 10% extra volume to avoid volume loss due to pipetting.
- g) After a few minutes of incubation, the instrument will pause and allow the user to add mixture prepared in step f. Open the cover and remove the plate from the instrument and add appropriate volume of Isopropanol (not provided) and Magnetic Beads C mixture that was prepared in step f to Column 1 and/or Column 7.

- h) Put the plate back into the instrument and press OK. The instrument will continue to perform the extraction.
- i) To rebind the sample to the magnetic beads, prepare a mixture of 200 µL **Buffer RL** and 300 µL Isopropanol for the number of samples extracted. (For example, if RNA is to be extracted from 10 samples, and mix 2 mL of **Buffer RL** with 3 mL **Isopropanol**). Make 10% extra volume to avoid volume loss due to pipetting.
- j) After a few minutes, the instrument will pause and allow the user to add mixture prepared in **step** i. Open the cover and remove the plate from the instrument and 500 μL of **Buffer RL** and **Isopropanol** mixture that was prepared in **step** i to **Column 2** and/or **Column 8**.
- k) Put the plate back into the instrument and press OK. The instrument will continue to perform the extraction.
- I) The instrument will prompt after the procedure is complete. Remove the elution plate from the instrument and the RNA is now ready for further downstream processing.
- m) Transfer the elution from **Column 6 and/or Column 12** and transfer it to a fresh plate for storage. The purified RNA sample may be stored at -20°C for a few days. It is recommended that samples be placed at -70°C for long-term storage.

2.5 Total RNA Purification - Manual method

- a) Vortex the lysate in RNase-free microcentrifuge tube (not provided) and rotate the tubes on a tube revolver rotator at 10 RPM at room temperature for 5 minutes. In case the tube revolver rotator is not available, incubate the tubes at room temperature for 5 minutes and occasionally vortex the tube.
- b) Add Isopropanol and Magnetic Bead C as per Table 5 to the lysate and vortex the tubes.
- c) Rotate the tubes on a tube revolver rotator at 10 RPM at room temperature for 5 minutes. In case the tube revolver rotator is not available, incubate the tubes at room temperature for 5 minutes and occasionally vortex the tube.
- d) Assemble a magnetic separation rack and place the sample tube in the magnetic rack. Allow to sit for 1 minute. Incubate for more time if the solution is not clear.
- e) Aspirate and discard supernatant without touching the magnetic beads.

Optional DNase Treatment

Norgen's Total RNA Purification Kit (Magnetic Bead System) isolates total RNA with minimal amounts of genomic DNA contamination. However, an optional **DNA Removal Protocol** for manual extraction is provided in **Appendix A** for maximum removal of residual DNA that may affect sensitive downstream applications. It is recommended that Norgen's RNase-Free DNase I Kit (Product # 25710) be used for this step. This step should be performed at this point in the protocol.

- a) Remove the tube from the magnetic rack and add 600 μ L of **Wash Solution A (ensure ethanol was added).**
- b) Rotate the tubes on a tube revolver rotator at 10 RPM at room temperature for 2 minutes. In case the tube revolver rotator is not available, incubate the tubes at room temperature for 2 minutes and gently vortex the tube, occasionally.
- c) Place the sample tube on the magnetic rack and allow to sit for 1 minute. Incubate for more time if the solution is not clear.
- d) Aspirate and discard supernatant without touching the magnetic beads.
- e) Repeat **steps 2f to 2i** two more times, for a total of three washes.
- f) Remove as much of the wash solution as possible and air-dry the beads at room temperature for 5 to 10 min. (Do not let the magnetic beads to over dry as this may affect the RNA yield).

- g) Remove the tube from the magnetic rack and add 50 100 μL of Elution Solution A.
- h) Mix by pipetting and incubate at 50°C for 2 minutes.
- i) Incubate the tubes at room temperature for 2 minutes and occasionally gently mixing by vortexing the tube. Briefly spin down the tubes to bring the contents at the bottom of the tube.
- j) Place the tube on the magnetic rack. Allow to sit for 1 minute. Incubate for more time if the solution is not clear.
- k) Carefully transfer the elution to a fresh 1.7 mL elution tube (provided) without touching the magnetic beads.
- I) The purified RNA sample may be stored at -20°C for a few days. It is recommended that samples be placed at -70°C for long-term storage.

Appendix A

Protocol for Optional DNA Removal

Norgen's Total RNA Purification Kit isolates total RNA with minimal amounts of genomic DNA contamination. However, an optional protocol is provided below for maximum removal of residual DNA that may affect sensitive downstream applications. It is recommended that Norgen's RNase- Free DNase I Kit (Product # 25710) be used for this step.

For every DNase reaction to be performed, prepare a DNase I Mix of 15 μL of DNase I and 100 μL of Enzyme Incubation Buffer A using Norgen's RNase-Free DNase I Kit (Product# 25710). Mix gently by inverting the tube a few times. DO NOT VORTEX.

Note: If using an alternative **DNase I**, prepare a working stock of 0.25 Kunitz unit/ μ L RNase-free DNase I solution according to the manufacturer's instructions. A 100 μ L aliquot is required for each sample to be treated.

- 2. Perform appropriate binding up to **step d of section 2.5** of the manual extraction and remove the tube from the magnetic rack.
- Add 100 μL of **DNase I Mix** that to the sample tube and mix by pipetting and incubate at 25-30°C for 15 minutes.
- 4. Add 200 µL **Buffer RL** and 300 µL **Isopropanol** to each sample.
- 5. Rotate the tubes on a tube revolver rotator at 10 RPM at room temperature for 5 minutes. In case the tube revolver rotator is not available, incubate the tubes at room temperature for 5 minutes and gently vortex the tube, occasionally.
- 6. Place tubes on the magnetic rack. Allow to sit for 1 minute. Allow to sit for 1 minute. Incubate for more time if the solution is not clear.
- 7. Aspirate and discard supernatant without touching the magnetic beads.
- 8. Proceed to **step f** of **section 2.5** to continue manual extraction.

Troubleshooting Guide

Problem	Possible Cause	Solution and Explanation	
Magnetic beads were accidentally pipetted up with the supernatant.	The pipette tip was placed too close to the magnetic beads while pipetting	Return the magnetic beads and the supernatant back into the sample tube. Mix well, and place the tube back onto the magnetic separation rack for the specified time. Carefully remove the supernatant without touching the magnetic beads.	
	Incomplete lysis of cells	Ensure that correct lysis protocol was applied to the sample.	
The yield of RNA	Amount of magnetic beads added was not sufficient	Ensure that the Magnetic Beads C is mixed well prior to use to avoid any inconsistency in RNA isolation.	
	RNA concentration in the cell or tissue sample being used is low.	Some samples contain very little target RNA. This varies from individual to individual based on numerous variables.	
Very gelatinous prior to adding the	The lysate solution mixture is not homogeneous	To ensure a homogeneous solution, vortex for 10-15 seconds before adding the magnetic beads to the lysate.	
Magnetic bead and Isopropanol	Maximum number of cells or amount of tissue exceeds kit specifications	Refer to specifications to determine if amount of starting material falls within kit specifications.	
RNA does not perform well in downstream	DNA was not washed with Wash Solution A	Traces of salt from the binding step may remain in the sample if the magnetic beads are not washed with Wash Solution A. Salt may interfere with downstream applications, and thus must be washed from the magnetic beads.	
applications.	Ethanol carryover	Ensure that the drying step after the last wash steps is performed, in order to remove traces of ethanol prior to elution. Ethanol is known to interfere with many downstream applications.	
DNA is present in eluted RNA.	DNA is coeluted with the RNA	Carry out a digestion with DNase I on the elution if the DNase present will interfere with downstream applications. Refer to manufacturer's instructions regarding amount of enzyme to use, optimal incubation time and temperature.	

Related Products (Required for Automated Workflow)	Product #
Adhesive Tape	28394

3430 Schmon Parkway, Thorold, ON Canada L2V 4Y6 Phone: (905) 227-8848 Fax: (905) 227-1061

Toll Free in North America: 1-866-667-4362

©2024 Norgen Biotek Corp.