

Plasma/Serum Nucleic Acid




Isolation Kits for 50 µL - 5 mL




- Isolate cell-free circulating and viral RNA & DNA
- No phenol extractions, no carrier RNA
- Concentrate into 50 µL to 100 µL
- High quality & purity RNA in 15-40 minutes
- Sensitive, linear and inhibitor-free nucleic acids
- For use with Norgen's cf-DNA/cf-RNA Preservative Tubes


Robust Plasma Volumes



Excellent Linearity



Inhibitor-Free RNA



Nucleic Acids isolated using Norgen's kit were used as a template in qPCR and qRT-PCR to assess relative amount of DNA and RNA linearity and level of inhibition at various volumes.

High Quality Nucleic Acid Enabling

Methylation Epigenetics

NGS

Fetal DNA Screening

Microarrays

Biomarker Discovery

Viral Diagnostics

Ordering information

Description	Size	Cat. Number
Plasma/Serum Cell-Free Circulating and Viral Nucleic Acid Purification Mini Kit	50 preps	56300
Plasma/Serum Cell-Free Circulating and Viral Nucleic Acid Purification Midi Kit	20 preps	56400
Plasma/Serum Cell-Free Circulating and Viral Nucleic Acid Purification Maxi Kit	10 preps	56500

