

RNA Clean-Up and Concentration Micro-Elute Kit




Rapid and efficient clean-up and concentration of total RNA, including microRNA, without phenol from small input volumes




- ✓ Clean-up and concentrate total RNA (including miRNA) in minutes
- ✓ Ideal for concentrating RNA purified from exosomes, plasma, serum, urine and other bodily fluids into small elution volumes.
- ✓ Purify all sizes of RNA, from large mRNA down to miRNA
- ✓ Clean-up RNA from contaminants including enzymes, primers, nucleotides and endotoxins
- ✓ Efficient RNA clean-up from enzymatic reactions, including labelling and DNase treatments
- ✓ Purified RNA is suitable for a variety of downstream applications, including Small RNA Sequencing


Excellent Concentration for RNA: Ideal for NGS Applications




High Recovery of Concentrated RNA




Excellent Quality of Concentrated RNA



Concentration of Total RNA



Concentration of miRNA



About this kit

Norgen's RNA Clean-Up and Concentration Micro-Elute Kit provides a rapid method for the purification, cleanup and concentration of RNA for NGS Library preparation and excellent processing of up to 45 µg of RNA isolated using different methods including phenol/guanidine-based protocols, and from various upstream enzymatic reactions such as DNase treatment, labeling and in vitro transcription.

Ordering information

Description	Size	Cat. Number
RNA Clean-Up and Concentration Micro-Elute Kit	50 preps	61000
RNA Clean-Up and Concentration Kit	50 preps	23600
RNA Clean-Up and Concentration Kit	100 preps	43200
RNA Clean-Up and Concentration 96-Well Kit	2 x 96-well plates	25100
Total RNA Purification Kit	50 preps	17200
Total RNA Purification Kit	100 preps	37500

v3.0

