

***Streptococcus uberis* TaqMan Probe/Primer and Control Set - 100 Reactions**
Product # TM30810

Component	Cat# TM30810 (100 Reactions)
Streptococcus uberis Primer & Probe Mix	1 x 280 µL
Streptococcus uberis Positive Control	1 x 150 µL
Nuclease-Free Water	1 x 1.25 mL

Applications

TaqMan Probe/Primer and Control Set for *Streptococcus uberis* detection using real-time PCR based on the use of TaqMan® technology.

Storage Conditions and Product Stability

Store all the components at -20°C upon arrival. All kit components can be stored for 2 years after the date of production without showing any reduction in performance.

Precautions and Disclaimers

This product is designed for research purposes only. It is not intended for human or diagnostic use.

Customer-Supplied Reagents and Equipment

- Real-Time PCR Instrument with multiple filter channel (FAM channel minimum)
- 2x PCR Master Mix
- **Optional** – Norgen's TaqMan 2X PCR Master Mix can be used for PCR validation purposes (Please see related products table below)

Procedures

Recommended Reaction Conditions:

- Streptococcus uberis Primer & Probe Mix (2 µL per 20 µL PCR reaction)
- Streptococcus uberis Positive Control (2-8 µL per 20 µL PCR reaction)

Description

The positive control is a cloned partial *Streptococcus uberis* DNA.

Recommended TaqMan PCR Assay

PCR Cycle	Step	Temperature	Duration
<i>Cycle 1</i>	Step 1	95°C	3 min
<i>Cycle 2 (40x)</i>	Step 1	95°C	15 sec
	Step 2	60°C	30 sec

Interpretation of Assay Results

* FAM (Target detection)	** HEX/VIC (PCR validation-optional)	Result
+	+	Positive
-	+	Negative
-	-	PCR inhibited

* FAM channel for target detection with a regular 2x PCR Master Mix

** **Optional** - HEX or VIC channel is only required when Norgen's TaqMan 2X PCR Master Mix (Cat.#28340) is used for PCR validation purpose.

Related Products	Product #
Streptococcus uberis TaqMan TaqMan PCR Kit	TM30850
Plasma/Serum cfc-DNA Advanced Purification Kit	68000
TaqMan 2X PCR Master Mix	28340

Technical Support

Contact our Technical Support Team between the hours of 8:30 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362. Technical support can also be obtained from our website or through email at techsupport@norgenbiotek.com.

Norgen Biotek Corp.
3430 Schmon Parkway, Thorold, ON Canada L2V 4Y6
Phone: (905) 227-8848
Fax: (905) 227-1061
Toll Free in North America: 1-866-667-4362