

Ad5 Type C TaqMan Probe/Primer and Control Set - 100 Reactions Product # TM42610

Component	Product # TM42610 (100 Reactions)
Ad5 Type C Primer & Probe Mix	1 x 280 μL
Ad5 Type C Positive Control	1 x 150 μL
Nuclease-Free Water	1 x 1.25 mL

Applications

TaqMan Probe/Primer and Control Set for Ad5 Type C detection using real-time PCR based on the use of TaqMan® technology.

Storage Conditions and Product Stability

Store all the components at -20°C upon arrival. All kit components can be stored for 2 years after the date of production without showing any reduction in performance.

Precautions and Disclaimers

This product is designed for research purposes only. It is not intended for human or diagnostic use.

Customer-Supplied Reagents and Equipment

- Real-Time PCR Instrument with multiple filter channel (FAM channel minimum)
- 2x PCR Master Mix
- Optional Norgen's TaqMan 2X PCR Master Mix can be used for PCR validation purposes (Please see related products table below)

Procedures

Recommended Reaction Conditions:

- Ad5 Type C Primer and Probe Mix (2 μL per 20 μL PCR reaction)
- Ad5 Type C Positive Control (2-8 μL per 20 μL PCR reaction)

Description

The positive control is a synthetic oligonucleotide containing partial Ad5 Type C.

Recommended TaqMan one step PCR Assay

One Step PCR Cycle	Step	Temperature	Duration
Cycle 1	Step 1	95°C	3 min
Cycle 2 (40x)	Step 1	95°C	15 sec
	Step 2	60°C	30 sec

Interpretation of Assay Results

* FAM (Target detection)	** HEX/VIC (optional for PCR validation)	Result
+	+	Positive
-	+	Negative
-	-	PCR inhibited

^{*} FAM channel for target detection with a regular 2x PCR Master Mix

^{**} **Optional** - HEX or VIC channel is only required when Norgen's TaqMan 2X PCR Master Mix (Cat.#28340) is used for PCR validation purpose.

Related Products	Product #
Ad5 Type C TaqMan PCR Kit	TM42650
TaqMan 2X PCR Master Mix	28340
2X PCR Master Mix	28007, 28008, 28009

Technical Support

Contact our Technical Support Team between the hours of 9:00 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362. Technical support can also be obtained from our website or through email at techsupport@norgenbiotek.com.

© 2024 Norgen Biotek Corp. PITM42610-1