

Toxoplasma gondii TaqMan Probe/Primer and Control Set - 100 Reactions Product # TM44710

Component	Product # TM44710 (100 Reactions)
Toxoplasma gondii Primer & Probe Mix	1 x 280 μL
Toxoplasma gondii Positive Control	1 x 150 μL
Nuclease-Free Water	1 x 1.25 mL

Applications

TaqMan Probe/Primer and Control Set for *Toxoplasma gondii* detection using real-time PCR based on the use of TaqMan® technology.

Storage Conditions and Product Stability

Store all the components at -20°C upon arrival. All kit components can be stored for 2 years after the date of production without showing any reduction in performance.

Precautions and Disclaimers

This product is designed for research purposes only. It is not intended for human or diagnostic use.

Customer-Supplied Reagents and Equipment

- Real-Time PCR Instrument with multiple filter channel (FAM channel minimum)
- 2x PCR Master Mix
- Optional Norgen's TaqMan 2X PCR Master Mix can be used for PCR validation purposes (Please see related products table below)

Procedures

Recommended Reaction Conditions:

- Toxoplasma gondii Primer Mix (2 μL per 20 μL PCR reaction)
- Toxoplasma gondii Positive Control (2-8 μL per 20 μL PCR reaction)

Description

The positive control is a cloned partial *Toxoplasma gondii* DNA.

Recommended TagMan PCR Assay

PCR Cycle	Step	Temperature	Duration
Cycle 1	Step 1	95°C	3 min
Overla 2 (40x)	Step 1	95°C	15 sec
Cycle 2 (40x)	Step 2	60°C	30 sec

Interpretation of Assay Results

* FAM (Target detection)	** HEX/VIC (PCR validation-optional)	Result
+	+	Positive
-	+	Negative
-	-	PCR inhibited

^{*} FAM channel for target detection with a regular 2x PCR Master Mix

^{**} **Optional** - HEX or VIC channel is only required when Norgen's TaqMan 2X PCR Master Mix (Cat.#28340) is used for PCR validation purpose.

Related Products	Product #
Toxoplasma gondii TaqMan PCR Kit	TM44750
TaqMan 2X PCR Master Mix	28340
2X PCR Master Mix	28007, 28008, 28009

Technical SupportContact our Technical Support Team between the hours of 8:30 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362. Technical support can also be obtained from our website or through email at techsupport@norgenbiotek.com.

© 2022 Norgen Biotek Corp. PITM44710-2