

Pospiviroid TaqMan Probe/Primer and Control Set - 100 Reactions

Product # TM53710

Component	Product # TM53710 (100 Reactions)
Pospiviroid Primer & Probe Mix	1 x 280 µL
Pospiviroid Positive Control	1 x 150 µL
Nuclease-Free Water	1 x 1.25 mL

Applications

TaqMan Probe/Primer and Control Set for Pospiviroid detection using real-time PCR based on the use of TaqMan® technology

Storage Conditions and Product Stability

Store all the components at -20°C upon arrival. All kit components can be stored for 1 year after the date of production without showing any reduction in performance.

Precautions and Disclaimers

This product is designed for research purposes only. It is not intended for human or diagnostic use.

Customer-Supplied Reagents and Equipment

Appropriate Real-Time PCR Instrument with FAM and HEX filter channel

Procedures

Recommended Reaction Conditions:

- Pospiviroid Primer Mix (2 µL per 20 µL PCR reaction)
- Pospiviroid Positive Control (2-8 µL per 20 µL PCR reaction)

Description

The positive control is a cloned partial Pospiviroid IVT.

Recommended TaqMan one step RT-PCR Assay

One Step RT-PCR Cycle	Step	Temperature	Duration
<i>Cycle 1</i>	Step 1	50°C	30 min
<i>Cycle 2</i>	Step 1	95°C	3 min
<i>Cycle 3 (40x)</i>	Step 1	95°C	15 sec
	Step 2	60°C	30 sec

Technical Support

Contact our Technical Support Team between the hours of 8:30 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362. Technical support can also be obtained from our website or through email at techsupport@norgenbiotek.com.